



Proteogenix

Proteogenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

Proteogenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD62L/SELL Monoclonal Antibody
PRODUCT REFERENCE	ARO-A15777
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4
TARGET	CD62L/SELL

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	CD62 antigen-like family member L, Leukocyte adhesion molecule 1, CD62L, LYAM1, TQ1, SELL, Leukocyte-endothelial cell adhesion molecule 1, L-selectin, LAM-1, gp90-MEL, LNHR, LECAM1, Lymph node homing receptor, Leukocyte surface antigen Leu-8

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY