

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD63 Monoclonal Antibody
PRODUCT REFERENCE	ARO-A15786
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	CD63

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Ocular melanoma-associated antigen, TSPAN30, CD63 antigen, LAMP-3, Tetraspanin-30, OMA81H, Lysosomal-associated membrane protein 3, Tspan-30, Granulophysin, MLA1, Melanoma-associated antigen ME491, CD63