

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD66c/CEACAM6 Monoclonal Antibody
PRODUCT REFERENCE	ARO-A15798
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	CD66c/CEACAM6

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Normal cross-reacting antigen, Carcinoembryonic antigen-related cell adhesion molecule 6, CEACAM6, CD66c, NCA, Non-specific crossreacting antigen