

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD41/ITGA2B & CD61/ITGB3 Monoclonal Antibody
PRODUCT REFERENCE	ARO-A15815
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	CD41/ITGA2B & CD61/ITGB3

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	CD41, ITGAB, GPIIb, ITGA2B, Platelet membrane glycoprotein IIb, GP2B, Integrin alpha-IIb, GPalpha IIb, GPIIIa, Integrin beta-3, CD61, ITGB3, Platelet membrane glycoprotein IIIa, GP3A