



Proteogenix

Proteogenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

Proteogenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD44 Monoclonal Antibody
PRODUCT REFERENCE	ARO-A15824
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	CD44

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Hyaluronate receptor, MDU2, HUTCH-I, MDU3, Phagocytic glycoprotein I, CDw44, Epican, Phagocytic glycoprotein 1, ECMR-III, CD44, Hermes antigen, PGP-1, PGP-I, GP90 lymphocyte homing/adhesion receptor, LHR, CD44 antigen, MIC4, Extracellular matrix receptor

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY