

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD51/ITGAV Monoclonal Antibody
PRODUCT REFERENCE	ARO-A15831
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2, Kappa
TARGET	CD51/ITGAV

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	MSK8, VNRA, Vitronectin receptor, Integrin alpha-V, Vitronectin receptor subunit alpha, CD51, ITGAV, VTNR