

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD309/KDR/VEGFR-2 (1121B) Monoclonal Antibody
PRODUCT REFERENCE	ARO-A15868
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	CD309/KDR/VEGFR-2 (1121B)

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	FLK-1, Kinase insert domain receptor, VEGFR2, Fetal liver kinase 1, CD309, Vascular endothelial growth factor receptor 2, Protein-tyrosine kinase receptor flk-1, FLK1, VEGFR-2, KDR