

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD335/NCR1/NKp46 Monoclonal Antibody
PRODUCT REFERENCE	ARO-A15881
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	CD335/NCR1/NKp46

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	LY94, Natural cytotoxicity triggering receptor 1, Natural killer cell p46-related protein, NK cell-activating receptor, NK-p46, CD335, NKp46, NCR1, hNKp46, Lymphocyte antigen 94 homolog