

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD365/KIM-1/TIM-1/HAVCR1 Monoclonal Antibody
PRODUCT REFERENCE	ARO-A15891
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	CD365/KIM-1/TIM-1/HAVCR1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	TIM, Kidney injury molecule 1, T-cell immunoglobulin mucin receptor 1, KIM1, KIM-1, TIM1, TIMD1, HAVCR1, T-cell immunoglobulin and mucin domain-containing protein 1, T-cell membrane protein 1, CD365, Hepatitis A virus cellular receptor 1, HAVcr-1, TIMD-1,