

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Human CD366/HAVCR2/TIM-3 (SAA0101) Monoclonal Antibody
<b>PRODUCT REFERENCE</b>	ARO-A15894
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG4, Kappa
<b>TARGET</b>	CD366/HAVCR2/TIM-3 (SAA0101)

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Homo sapiens (Human)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	PBS buffer (pH 7.4)
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	T-cell immunoglobulin mucin receptor 3, HAVcr-2, TIMD-3, T-cell membrane protein 3, HAVCR2, Hepatitis A virus cellular receptor 2, CD366, TIMD3, TIM-3, TIM3, T-cell immunoglobulin and mucin domain-containing protein 3