

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD274/PD-L1/B7-H1 (22C3) Monoclonal Antibody
PRODUCT REFERENCE	ARO-A15933
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	CD274/PD-L1/B7-H1 (22C3)

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	B7-H1, Programmed cell death 1 ligand 1, PDCD1 ligand 1, PDCD1L1, B7 homolog 1, PDCD1LG1, PDL1, hPD-L1, Programmed death ligand 1, B7H1, PD-L1, CD274