

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD218b/IL18RAP Monoclonal Antibody
PRODUCT REFERENCE	ARO-A15979
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	CD218b/IL18RAP

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	AcPL, Accessory protein-like, Interleukin-1 receptor 7, IL-18R-beta, CD218 antigen-like family member B, IL-1RAcPL, Interleukin-18 receptor accessory protein, CDw218b, IL-1R-7, IL1R7, IL-1R7, IL-18RAcP, CD218b, IL-1R accessory protein-like, IL18RAP, Inter