

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD202b/TEK Monoclonal Antibody
PRODUCT REFERENCE	ARO-A15991
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	CD202b/TEK

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	VMCM, TIE2, Tunica interna endothelial cell kinase, CD202b, Tyrosine-protein kinase receptor TIE-2, TEK, Angiopoietin-1 receptor, VMCM1, Endothelial tyrosine kinase, Tyrosine-protein kinase receptor TEK, p140 TEK, Tyrosine kinase with Ig and EGF homology