

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD158a/b1/b2 Monoclonal Antibody
PRODUCT REFERENCE	ARO-A16002
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4, Kappa
TARGET	CD158a/b1/b2

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	CD158 antigen-like family member A, KIR2DL1, Killer cell immunoglobulin-like receptor 2DL1, Natural killer-associated transcript 1, NKAT-1, MHC class I NK cell receptor, CD158A, p58 natural killer cell receptor clones CL-42/47.11, CD158a, p58.1 MHC class-