

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD161/KLRB1/NK1.1 Monoclonal Antibody
PRODUCT REFERENCE	ARO-A16014
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	CD161/KLRB1/NK1.1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	CLEC5B, KLRB1, CD161, HNKRP-1a, C-type lectin domain family 5 member B, NKRP1A, Natural killer cell surface protein P1A, Killer cell lectin-like receptor subfamily B member 1, NKR-P1A