

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD182/CXCR2 Monoclonal Antibody
PRODUCT REFERENCE	ARO-A16015
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	CD182/CXCR2

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	IL8RB, GRO/MGSA receptor, C-X-C chemokine receptor type 2, CD182, CXCR-2, IL-8R B, CXC-R2, IL-8 receptor type 2, High affinity interleukin-8 receptor B, CDw128b, CXCR2