

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Human CD184/CXCR4 Monoclonal Antibody
<b>PRODUCT REFERENCE</b>	ARO-A16024
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1
<b>TARGET</b>	CD184/CXCR4

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Homo sapiens (Human)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	PBS buffer (pH 7.4)
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	Stromal cell-derived factor 1 receptor, LCR1, HM89, LPS-associated protein 3, LAP-3, CXC-R4, CXCR-4, SDF-1 receptor, CD184, NPYRL, LESTR, Fusin, FB22, C-X-C chemokine receptor type 4, CXCR4, Lipopolysaccharide-associated protein 3, Leukocyte-derived seven