

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD167a/DDR1 Monoclonal Antibody
PRODUCT REFERENCE	ARO-A16026
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	CD167a/DDR1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Protein-tyrosine kinase RTK-6, Protein-tyrosine kinase 3A, NTRK4, EDDR1, CAK, Tyrosine kinase DDR, Mammary carcinoma kinase 10, TRKE, Epithelial discoidin domain-containing receptor 1, CD167a, TRK E, Epithelial discoidin domain receptor 1, HGK2, CD167 ant