

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD16b/FCGR3B Monoclonal Antibody
PRODUCT REFERENCE	ARO-A16032
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	CD16b/FCGR3B

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	IGFR3, FCGR3B, CD16b, Low affinity immunoglobulin gamma Fc region receptor III-B, Fc-gamma RIII, IgG Fc receptor III-1, FCGR3, CD16B, Fc-gamma RIII-beta, Fc-gamma RIIIb, FCG3, FcRIII, FcR-10, FcRIIIb