

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|--------------------------------------|
| PRODUCT NAME | Human CD191/CCR1 Monoclonal Antibody |
| PRODUCT REFERENCE | ARO-A16041 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG1, kappa |
| TARGET | CD191/CCR1 |

PRODUCT INFORMATION

| | |
|--------------------------------|---|
| ORIGIN SPECIES | <i>Homo sapiens (Human)</i> |
| FORM | Liquid |
| BUFFER | PBS buffer (pH 7.4) |
| STABILITY & STORAGE | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| ALIAS | CMKR1, SCYAR1, Macrophage inflammatory protein 1-alpha receptor, CMKBR1, MIP-1alpha-R, HM145, CC-CKR-1, RANTES-R, LD78 receptor, CCR-1, C-C chemokine receptor type 1, CCR1, CD191, C-C CKR-1 |