

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD175/TNC Monoclonal Antibody
PRODUCT REFERENCE	ARO-A16044
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2b
TARGET	CD175/TNC

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Myotendinous antigen, Neuronectin, GMEM, Tenascin, Tenascin-C, TN-C, TNC, TN, Cytotactin, GP 150-225, Glioma-associated-extracellular matrix antigen, HXB, Hexabrachion, JI