

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD140a/PDGFR α (SAA0056) Monoclonal Antibody
PRODUCT REFERENCE	ARO-A16045
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	CD140a/PDGFR α (SAA0056)

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	PDGFR2, Platelet-derived growth factor receptor alpha, RHEPDGFR α , Alpha platelet-derived growth factor receptor, PDGFR-2, Platelet-derived growth factor alpha receptor, Platelet-derived growth factor receptor 2, PDGF-R-alpha, Alpha-type platelet-derived g