

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Human CD178/FASLG/TNFSF6 (NOK1/NOK-1) Monoclonal Antibody
<b>PRODUCT REFERENCE</b>	ARO-A16047
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG2a, kappa
<b>TARGET</b>	CD178/FASLG/TNFSF6 (NOK1/NOK-1)

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Homo sapiens (Human)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	PBS buffer (pH 7.4)
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	Receptor-binding FasL ectodomain, APL, APT1LG1, FasL ICD, Fas antigen ligand, CD178, sFasL, CD95L, APTL, TNFSF6, Apoptosis antigen ligand, SPA, FASLG, FasL, SPPL2A-processed FasL form, Tumor necrosis factor ligand superfamily member 6, Soluble Fas ligand,