

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human ERBB4/Her4 Monoclonal Antibody
PRODUCT REFERENCE	ARO-A16051
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a
TARGET	ERBB4/Her4

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	HER4, Tyrosine kinase-type cell surface receptor HER4, p180erbB4, s80HER4, ERBB4, 4ICD, Receptor tyrosine-protein kinase erbB-4, Proto-oncogene-like protein c-ErbB-4, E4ICD