

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human ERBB3/HER3 (SAA0123) Monoclonal Antibody
PRODUCT REFERENCE	ARO-A16054
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2-lambda
TARGET	ERBB3/HER3 (SAA0123)

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	ERBB3, HER3, Proto-oncogene-like protein c-ErbB-3, Tyrosine kinase-type cell surface receptor HER3, Receptor tyrosine-protein kinase erbB-3