

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD134/TNFRSF4/OX40 (SAA0051) Monoclonal Antibody
PRODUCT REFERENCE	ARO-A16062
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2, Kappa
TARGET	CD134/TNFRSF4/OX40 (SAA0051)

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	ACT35 antigen, TXGP1L, TAX transcriptionally-activated glycoprotein 1 receptor, Tumor necrosis factor receptor superfamily member 4, TNFRSF4, OX40L receptor, CD134