

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human EPHA3 Monoclonal Antibody
PRODUCT REFERENCE	ARO-A16063
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	EPHA3

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	HEK, Eph-like tyrosine kinase 1, EK4, EPHA3, Tyrosine-protein kinase receptor ETK1, Tyrosine-protein kinase TYRO4, Ephrin type-A receptor 3, hEK4, ETK, TYRO4, ETK1, Human embryo kinase, EPH-like kinase 4