

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD120b/TNFRSF1B/TNFR2 Monoclonal Antibody
PRODUCT REFERENCE	ARO-A16067
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	CD120b/TNFRSF1B/TNFR2

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	TNF-R2, TNFR-II, p80 TNF-alpha receptor, Etanercept, TNFBR, TBP-2, TBPII, Tumor necrosis factor receptor superfamily member 1B, CD120b, p75, TNFRSF1B, Tumor necrosis factor receptor 2, TNF-RII, TNFR2, Tumor necrosis factor receptor type II