

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD121a/IL1R1 Monoclonal Antibody
PRODUCT REFERENCE	ARO-A16070
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	CD121a/IL1R1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	IL-1RT1, IL-1R-alpha, mL-1R1, CD121a, IL-1R-1, sIL-1RI, IL1RT1, IL1R1, sIL-1R1, IL1R, CD121 antigen-like family member A, Interleukin-1 receptor type I, Interleukin-1 receptor alpha, mL-1RI, IL1RA, IL-1RT-1, Interleukin-1 receptor type 1, p80