

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD136/MST1R Monoclonal Antibody
PRODUCT REFERENCE	ARO-A16071
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	CD136/MST1R

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	CD136, Macrophage-stimulating protein receptor, PTK8, MST1R, p185-Ron, RON, CDw136, MSP receptor, Protein-tyrosine kinase 8