

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Human CD124/IL4R Monoclonal Antibody
<b>PRODUCT REFERENCE</b>	ARO-A16081
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG4, Kappa
<b>TARGET</b>	CD124/IL4R

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Homo sapiens (Human)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	PBS buffer (pH 7.4)
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	IL-4RA, IL-4-binding protein, IL-4 receptor subunit alpha, Soluble IL-4R-alpha, IL4-BP, IL4R, IL-4R subunit alpha, IL-4R-alpha, CD124, Interleukin-4 receptor subunit alpha, IL4RA, Soluble IL-4 receptor subunit alpha, sIL4Ralpha/prot