

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human FOLR1 (SAA0118) Monoclonal Antibody
PRODUCT REFERENCE	ARO-A16083
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	FOLR1 (SAA0118)

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Folate receptor alpha, FOLR, FOLR1, Adult folate-binding protein, Folate receptor 1, FR-alpha, KB cells FBP, FBP, Ovarian tumor-associated antigen MOv18, Folate receptor, adult