

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD126/IL6R/IL-6RA (SAA0048) Monoclonal Antibody
PRODUCT REFERENCE	ARO-A16089
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	CD126/IL6R/IL-6RA (SAA0048)

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	sIL6R, Interleukin-6 receptor subunit alpha, CD126, IL-6R 1, IL-6R subunit alpha, Membrane glycoprotein 80, gp80, IL6R, IL-6RA, IL-6 receptor subunit alpha, IL-6R-alpha