

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Human CD126/IL6R/IL-6RA Monoclonal Antibody
<b>PRODUCT REFERENCE</b>	ARO-A16092
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG2, Kappa
<b>TARGET</b>	CD126/IL6R/IL-6RA

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Homo sapiens (Human)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	PBS buffer (pH 7.4)
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	sIL6R, Interleukin-6 receptor subunit alpha, CD126, IL-6R 1, IL-6R subunit alpha, Membrane glycoprotein 80, gp80, IL6R, IL-6RA, IL-6 receptor subunit alpha, IL-6R-alpha