

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CLEC5A/MDL1 Monoclonal Antibody
PRODUCT REFERENCE	ARO-A16094
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	CLEC5A/MDL1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	C-type lectin superfamily member 5, C-type lectin domain family 5 member A, CLEC5A, MDL-1, CLECSF5, MDL1, Myeloid DAP12-associating lectin 1