

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CLEC2D/OCIL/LLT1 Monoclonal Antibody
PRODUCT REFERENCE	ARO-A16097
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	CLEC2D/OCIL/LLT1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	CLAX, C-type lectin domain family 2 member D, Osteoclast inhibitory lectin, CLEC2D, Lectin-like transcript 1, Lectin-like NK cell receptor, OCIL, LLT-1, LLT1