



ProteoGenix

ProteoGenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

ProteoGenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|--------------------------------------|
| PRODUCT NAME | Human GOLM1/GP73 Monoclonal Antibody |
| PRODUCT REFERENCE | ARO-A16114 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG1 |
| TARGET | GOLM1/GP73 |

PRODUCT INFORMATION

| | |
|--------------------------------|---|
| ORIGIN SPECIES | <i>Homo sapiens (Human)</i> |
| FORM | Liquid |
| BUFFER | PBS buffer (pH 7.4) |
| STABILITY & STORAGE | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| ALIAS | Golgi phosphoprotein 2, C9orf155, Golgi membrane protein GP73, Golgi membrane protein 1, GOLPH2, GOLM1 |

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY