

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human GZMB (scFv17) Monoclonal Antibody
PRODUCT REFERENCE	ARO-A16148
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	GZMB (scFv17)

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Cathepsin G-like 1, Granzyme-2, GZMB, Fragmentin-2, Human lymphocyte protein, HLP, CTLA-1, CGL1, CTSSL1, CSPB, Lymphocyte protease, T-cell serine protease 1-3E, Granzyme B, CTLA1, GRB, C11, Cytotoxic T-lymphocyte proteinase 2, SECT