

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human VTCN1/B7-H4 Monoclonal Antibody
PRODUCT REFERENCE	ARO-A16160
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	VTCN1/B7-H4

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	V-set domain-containing T-cell activation inhibitor 1, B7h.5, B7 homolog 4, B7-H4, T-cell costimulatory molecule B7x, B7H4, VTCN1, Protein B7S1, Immune costimulatory protein B7-H4