

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human TACSTD2/TROP2 (TP15-4) Monoclonal Antibody
PRODUCT REFERENCE	ARO-A16169
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	TACSTD2/TROP2 (TP15-4)

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Pancreatic carcinoma marker protein GA733-1, Tumor-associated calcium signal transducer 2, M1S1, TROP2, TACSTD2, Cell surface glycoprotein Trop-2, GA733-1, Membrane component chromosome 1 surface marker 1