

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|---|
| PRODUCT NAME | Human SORT1/Gp95/Sortilin Monoclonal Antibody |
| PRODUCT REFERENCE | ARO-A16178 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG1, kappa(silence) |
| TARGET | SORT1/Gp95/Sortilin |

PRODUCT INFORMATION

| | |
|--------------------------------|---|
| ORIGIN SPECIES | <i>Homo sapiens (Human)</i> |
| FORM | Liquid |
| BUFFER | PBS buffer (pH 7.4) |
| STABILITY & STORAGE | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| ALIAS | Neurotensin receptor 3, Gp95, Glycoprotein 95, Sortilin, 100 kDa NT receptor, NTR3, SORT1, NT3 |