

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Human CD116/CSF2RA Monoclonal Antibody
PRODUCT REFERENCE	ARO-A16200
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4, lambda
TARGET	CD116/CSF2RA

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	PBS buffer (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	CD116, CSF2R, GMR-alpha, CSF2RA, Granulocyte-macrophage colony-stimulating factor receptor subunit alpha, CSF2RY, GMCSFR-alpha, GM-CSF-R-alpha, CDw116