

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Pityrosporum furfur Mal f 1/Allergen Pit o 1 NANOBODY [®] , protein (SAA1212)
PRODUCT REFERENCE	ARO-A16216
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	Mal f 1/Allergen Pit o 1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Malassezia furfur</i> (Pityriasis versicolor infection agent) (<i>Pityrosporum furfur</i>)
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-Major allergen Mal f 1, anti-Allergen Pit o 1, anti-Mal f 1