

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-Vibrio cholerae serotype O1 higB-2 NANOBODY&reg, protein (SAA1029)
<b>PRODUCT REFERENCE</b>	ARO-A16226
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	VHH-8His-Cys-tag
<b>TARGET</b>	higB-2

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Vibrio cholerae serotype O1 (strain ATCC 39315 / El Tor Inaba N16961)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4)
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	anti-Toxin HigB-2, anti-higB-2, anti-VC_A0468