

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Castor bean Ricin NANOBODY [®] , protein (SAA0984)
PRODUCT REFERENCE	ARO-A16227
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	Ricin

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Ricinus communis (Castor bean)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-Ricin, anti-Ricin A chain, anti-3.2.2.22, anti-rRNA N-glycosidase, anti-Ricin B chain