

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Clostridium botulinum BoNT/B NANOBODY [®] , protein (SAA0939)
PRODUCT REFERENCE	ARO-A16246
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	BoNT/B

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Clostridium botulinum</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-Botulinum neurotoxin type B, anti-BoNT/B, anti-Bontoxilysin-B, anti-Botulinum neurotoxin B light chain, anti-LC, anti-3.4.24.69, anti-Botulinum neurotoxin B heavy chain, anti-HC, anti-botB