

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Vibrio cholerae serotype O1 ctxB/Cholera Toxin Subunit B NANOBODY [®] protein (SAA1345)
PRODUCT REFERENCE	ARO-A16257
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	ctxB/Cholera Toxin Subunit B

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Vibrio cholerae serotype O1 (strain ATCC 39315 / El Tor Inaba N16961)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-Cholera enterotoxin subunit B, anti-Cholera enterotoxin B chain, anti-Cholera enterotoxin gamma chain, anti-Cholera toxin, anti-CTxB, anti-toxB, anti-VC_1456