

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-BC2-tag (PDRKAAVSHWQQ) NANOBODY®, protein (SAA1411)
PRODUCT REFERENCE	ARO-A16260
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-His
TARGET	CTNNB1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-CTNNB, anti-Catenin beta-1, anti-CTNNB1, anti-Beta-catenin