

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CD369/CLEC7A NANOBODY®, protein (SAA1364)
PRODUCT REFERENCE	ARO-A16268
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	CD369/CLEC7A

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-Beta-glucan receptor, anti-CD369, anti-CLEC7A, anti-DECTIN1, anti-BGR, anti-C-type lectin superfamily member 12, anti-DC-associated C-type lectin 1, anti-CLECSF12, anti-Dectin-1, anti-C-type lectin domain family 7 member A, anti-Dendritic cell-associ