

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Vesicular stomatitis Indiana virus/VSIV NP/Nucleoprotein NANOBODY [®] reg, protein (SAA1125)
PRODUCT REFERENCE	ARO-A16281
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	NP/Nucleoprotein

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Vesicular stomatitis Indiana virus (strain San Juan) (VSIV)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4)
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	anti-Nucleoprotein, anti-NP, anti-Nucleocapsid protein, anti-Protein N, anti-N, anti-Vesicular stomatitis Indiana virus (strain San Juan) (VSIV)